
femBloc[®]

Blended Polymer

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
Instructions for Use **Mode d'emploi** **Instrucciones de uso**

English	2
Français	4
Español	6

Blended Polymer

FBP-200

INSTRUCTIONS FOR USE

1. DEVICE DESCRIPTION

FemBloc Blended Polymer is a synthetic pre-mixed blend of cyanoacrylate-based tissue adhesive designed for use with the FemBloc Delivery System.

2. INTENDED PURPOSE

FemBloc Permanent Birth Control is intended for non-surgical permanent birth control by occlusion of the fallopian tubes. The FemBloc Blended Polymer is intended to stimulate wound healing response, gradually degrade and exit the body, resulting in complete tissue in-growth in a small section of the fallopian tubes.

2.1 Intended Patient Population

FemBloc Permanent Birth Control is intended for use in women aged 21-50.

2.2 Indication for Use

FemBloc Permanent Birth Control is indicated for non-surgical permanent birth control by occlusion of the fallopian tubes.

2.3 Contraindications

- Uncertainty about the desire to end fertility.
- Known or suspected pregnancy; or at risk for pregnancy from unprotected intercourse earlier in current cycle.
- Any condition or medical treatment (e.g. systemic corticosteroids or chemotherapy) that compromises immune system.
- Known hypersensitivity to cyanoacrylate or formaldehyde.

2.4 Clinical Benefits*

The primary clinical benefit of the FemBloc Permanent Birth Control is bilateral occlusion of the fallopian tubes for permanent contraception.

Table 1 and Table 2 provide data from the Femasys sponsored clinical trials for FemBloc procedure through 5 years. Pregnancy rate is shown in Table 1 and data demonstrating the safety events are shown in Table 2.

* Liu, J. H., et al. (2025). FemBloc Non-Surgical Permanent Contraception for Occlusion of the Fallopian Tubes. J Gynecol Reprod Med, 9(1), 01-19.

Table 1. Pregnancy Rate Reported in Clinical Trials with FemBloc

Statistic	Result (N=51) ^a
Overall, % (n/N)	0% (0/51)
Exact 95% upper confidence bound for p ₁	0.057
One-sided p-value ^b	0.0426

^a N=229 total patients who underwent the FemBloc procedure, N=101 were in the Cohort of Interest, and N=51 patients were determined bilaterally occluded by the investigator. ^b P-values <0.05 were considered indicative of significant difference from 6% performance goal based on reported pregnancy rate of 5.7% for female surgical sterilization. (Garipey et.al, Fertility and Sterility, 2022)

Table 2. Adverse Events >5% Reported in Clinical Trials with FemBloc

Adverse Event (AE)	% (N=229) ^a
Serious AE (SAE)	0
Non-Serious AE, device- or procedure-related	
Spotting vaginal or uterine/vaginal bleeding	58.5
Pelvic or abdominal pain/uterine cramps	55.9
Time of AE reporting from procedure	
1 day	61.1
2 days	8.5
3 days	16.2
>4 days	13.7
Non-Serious AE, device- or procedure-related	% (N=101) ^b
Bilateral/unilateral tubal patency	49.5

^a N=229 total patients who underwent the FemBloc procedure.

^b N=101 total patients who met trial eligibility, including Confirmation Test result 90 days after FemBloc procedure and N=51 patients who were determined bilaterally occluded by the investigator.

2.5 Intended Users

FemBloc Permanent Birth Control is intended to be used by licensed practitioners who are familiar with intrauterine transcervical techniques.

2.6 Performance Characteristics

The FemBloc Delivery System facilitates the accurate delivery of the FemBloc Blended Polymer into each uterine cavity cornu and fallopian tubes. Upon delivery and contact with biological tissue and tubal fluid, the Blended Polymer adhesive begins to polymerize. Through polymer formation, a pliable porous structure is formed that is adherent to the epithelial lining in a confined targeted area. The Blended Polymer induces irritation and, through a wound healing response, forms nonfunctional scar tissue while degrading and exiting the body through the vagina.

2.7 Device Lifetime

The FemBloc Blended Polymer has been shown to maintain performance through the procedure and stated shelf life. The majority of the delivered Blended Polymer is in the uterine cavity cornu, lined with endometrium, facilitating degradation during menses. The surface bonds of the Blended Polymer adhesive are constantly exposed to moisture and secretions, further aiding degradation. Excretion from the uterine cavity follows the direction of tubal fluid flow and pro-uterine ciliary beat, ultimately exiting through the vagina. Three months post the FemBloc procedure, a Confirmation Test is performed to assess for bilateral tubal occlusion with saline-air contrast that flushes the uterine cavity for a functional device lifetime of three months.

NOTE: MRI Safe as Blended Polymer poses no known hazards when used within any MRI environment and contains no metal.

INSTRUCTIONS FOR USE

Ensure use of aseptic technique throughout procedure.

3. WARNINGS & PRECAUTIONS

- FemBloc should be performed during: 1) follicular phase of the menstrual cycle (not to exceed day 11 from first day of menses) after bleeding has stopped, or 2) any time she is not bleeding and protected from pregnancy if on hormonal contraceptive.
- FemBloc should not be placed if there is active uterine bleeding, fluid in uterine cavity, or blood in vagina as bleeding/fluid may affect Blended Polymer delivery.
- Alternative contraception must be used until the final FemBloc Confirmation Test is performed. If Confirmation test does not document occlusion or is inconclusive, patient must use an alternative form of birth control.
- FemBloc should be considered irreversible.
- Sterile device. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize.
- Do not use after the use-by date.

4. CLINICAL INFORMATION

The Summary of Safety and Clinical Performance (SSCP) can be found at <http://ec.europa.eu/tools/eudamed>.

5. POTENTIAL ADVERSE EVENTS

- Hypersensitivity or allergic reactions to cyanoacrylates, such as urticaria, shortness of breath and anaphylactic shock.

6. PATIENT COUNSELING INFORMATION

Licensed practitioners should review the instructions for use for applicable information to be shared with the patient as an implant card will not be provided as the Blended Polymer will degrade and leave the body. The patient should be counseled to use alternative contraception (except IUD) until the Confirmation Test confirms bilateral tubal occlusion.

7. STORAGE

The FemBloc Blended Polymer should be stored in its carton at room temperature and humidity. The Blended Polymer pouch should only be opened when intending to use.

8. DEVICE MATERIALS

The FemBloc Blended Polymer is provided in two containers, each containing 0.9 mL of pre-mixed methoxy-n-propyl and methoxy-iso-propyl cyanoacrylate-based adhesive within the cartridge that is inserted into the FemBloc Delivery System for delivery of 0.6 mL of Blended Polymer directed towards each fallopian tube.

FBP-200 • Blended Polymer

FemBloc® Permanent Birth Control

9. INSTRUCTIONS FOR USE

9.1 Required Accessories

- FemBloc Delivery System

9.2 Device Preparation

- Visually inspect the sterile barrier system (SBS) for breaches of package integrity prior to aseptic presentation.
- Remove Blended Polymer cartridge from foil package by tearing open pre-cut slit.
- Remove cap from cartridge (Figure 1).

NOTE: If cap is not attached to cartridge, **Do Not Use** Blended Polymer. Open a **NEW** Blended Polymer package.

- Insert Blended Polymer cartridge into FemBloc Delivery System cartridge compartment by inserting barrel tips into holes in the compartment at an angle and then pressing down until cartridge snaps into place (Figure 2).

NOTE: Once cartridge is inserted, it becomes permanently contained within the FemBloc Delivery System and cannot be removed.


Figure 1


Figure 2

9.3 FemBloc Procedure

- Follow FemBloc Delivery System Instructions for Use (IFU) for FemBloc procedure preparation, Delivery System placement, and procedure completion.
- Use transabdominal ultrasound to assess for proper FemBloc Delivery System placement prior to Blended Polymer delivery.

9.4 Post Procedure Instructions and Patient Follow-Up

- Ensure patient uses alternative contraception (except IUD) until the final Confirmation Test (M-Sono HSG). Counsel patient on risk of pregnancy due to non-compliance.
- Schedule patient for Confirmation Test at least three (3) months after the FemBloc procedure.

10. DISPOSAL

Dispose of FemBloc Permanent Birth Control in accordance with applicable laws, regulations, and hospital procedures regarding biohazardous medical devices.

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
Tel: (31) (0) 70 345-8570

VIGILANCE REPORTING

Serious incidents occurring in the European Union in relation to FemBloc Blended Polymer should be reported to the Authorized Representative and/or the Manufacturer as indicated below:

Authorized Representative: Emergo – EmergoVigilance@ul.com

Manufacturer: Femsys Inc. – Vigilance@femsys.com

User/Patient: Serious incidents should be reported to the manufacturer and corresponding competent authority of the Member State.


Manufacturer
Femsys Inc.
3950 Johns Creek Court, Suite 100
Suwanee, GA 30024 USA
Tel: (877) 336-2562
Fax: (404) 581-5903
www.femsys.com
Patent Listing: www.femsys.com/patents

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Fabricante: Femsys Inc. – Vigilance@femasys.com
Femasys Inc.