
femBloc®

Delivery System

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Sistema de implantación


Instructions for Use **Mode d'emploi** **Instrucciones de uso**

English	3
Français	7
Español	11

Delivery System

Instructions for Use

DEVICE DESCRIPTION

FemBloc Delivery System consists of two balloon catheters (not made with natural rubber latex), insertion tube, and device handle.

INDICATIONS FOR USE

FemBloc Delivery System is intended to deliver FemBloc Blended Polymer into each uterine cavity cornu and both fallopian tubes.

The clinical indication: non-surgical permanent birth control by occlusion of the fallopian tubes.

STORAGE

Store in a cool, dry place.

HOW SUPPLIED

- Sterile for single use only.
- Sterilized with Ethylene Oxide.

REQUIRED ACCESSORY

- FemBloc Permanent Birth Control Blended Polymer (cartridge)

RECOMMENDED ANCILLARY STERILE SUPPLIES

- Vaginal speculum
- Tenaculum (straight)
- Sound

CONTRAINDICATIONS

- Uncertainty about the desire to end fertility.
- Known or suspected pregnancy; or at risk for pregnancy from unprotected intercourse earlier in current cycle.
- Known endometrial or myometrial conditions (e.g. submucous leiomyoma), uterine anomaly (e.g. unicornuate, bicornuate, arcuate, septate, or didelphic), or uterine position (e.g. retroflexion or antelexion) that would interfere with insertion tube midline fundal placement, access to uterine cornual region, or lateral deployment of the catheters.
- Any condition which may prohibit proper visualization of the cervix or not allow the uterus to be appropriately instrumented.
- Active upper or lower genital tract infection.
- Less than 6 weeks post partum/post pregnancy termination at time of FemBloc procedure.

WARNINGS & PRECAUTIONS

- FemBloc should be performed during: 1) follicular phase of the menstrual cycle (not to exceed day 11 from first day of menses) after bleeding has stopped, or 2) any time she is not bleeding and protected from pregnancy if on hormonal contraceptive.
- FemBloc should not be placed if there is active uterine bleeding, fluid in uterine cavity, or blood in vagina as bleeding/fluid may affect placement.
- Alternative contraception must be used until the final FemBloc Confirmation Test is performed. If Confirmation Test does not document occlusion or is inconclusive, patient must use an alternative form of birth control.
- FemBloc should be considered irreversible.
- Sterile device. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize.
- Dispose in accordance with all local, state, and national Medical/Hazardous Waste practices.

ADVERSE EVENTS

- Some patients may experience bleeding/spotting and/or cramps after the procedure. Most resolve the day of the procedure.

INSTRUCTIONS FOR USE

Ensure use of aseptic technique throughout procedure.

FDS-200 • Delivery System

FemBloc® Permanent Birth Control

I. DEVICE PREPARATION

Visually inspect the sterile barrier system (SBS) for breaches of package integrity prior to aseptic presentation.

1. PREPARE DELIVERY SYSTEM

- Place Delivery System package flat on a horizontal surface with white Tyvek side up and peel Tyvek from tray.
- Perform a function test of the balloons by following the steps in Figure 1. Check balloons by inflating for 10 seconds and then deflating.

Note: Balloons do not lay flat in fully deflated state; appear slightly inflated.

If either balloon does not function, **Do Not Use** device. Open a **NEW** FemBloc Delivery System package.


Figure 1

Balloon Test

Step A: Push plunger forward until resistance is felt, partially inflating balloons.

Hold for 10 Seconds.

Step B: SLOWLY retract plunger to deflate both catheter balloons.

2. REFER TO FEMBLOC BLENDED POLYMER INSTRUCTIONS FOR USE

- Insert Blended Polymer cartridge into Delivery System cartridge compartment by inserting barrel tips into holes in the compartment at an angle and then pressing down until cartridge snaps into place (Figure 2).

Figure 2



Figure 3

II. PROCEDURE PREPARATION

- Evaluate uterine cavity for presence of fluid/blood by inserting transvaginal probe and obtaining sagittal view.

Note: FemBloc Delivery System should not be placed if fluid/blood is detected in the uterine cavity.

- Grasp lip of cervix with tenaculum and control for bleeding prior to device placement.
- Measure depth of uterine cavity with sound.

Note: Minimize repeated tapping of the fundus for patient comfort.

- Set flange on insertion tube to correspond with uterine depth determined by sound (Figure 3).

III. DEVICE PLACEMENT

Device should be operated at all times with device top side up.


1. INSERT DELIVERY SYSTEM

- Prior to insertion, move slider back until catheters are fully retracted and contained within insertion tube (Figure 4).

Note: Do not leave catheters retracted within insertion tube for extended time as catheters may lose their shape (straighten).

- Securely hold device handle in one palm and orient in horizontal or uterine plane with uterine depth scale visible.

Figure 4

INSTRUCTIONS FOR USE

FDS-200 • Delivery System

FemBloc® Permanent Birth Control

WARNING: To reduce risk of uterine perforation, do not use force to overcome resistance if encountered during insertion.

- c. While maintaining gentle traction on cervix with tenaculum, advance insertion tube gently through cervix into uterine cavity. Stop advancing when fundal resistance is felt and confirm flange is touching external cervical os.

Note: If resistance is encountered, use dilation (<5 mm recommended) to overcome.

- d. The tip of insertion tube should now be at the desired fundal position for the remainder of the procedure (**Figure 5**). Fundal positioning is required for proper device placement.
- e. Secure tenaculum to device handle by fastening tenaculum stabilizer to the hook (**Figure 6**).

Note: Keep device handle oriented in horizontal or uterine plane AND maintain fundal position of insertion tube by confirming flange is positioned at external cervical os throughout entire procedure, as depicted in **Figure 5**.

2. DELIVER BLENDED POLYMER

- a. Confirm flange is located at external cervical os and continue to securely hold device handle in one palm throughout Blended Polymer delivery.
- b. With free hand, move slider completely forward until it clicks indicating the catheters are fully advanced (**Figure 8**).
- c. Confirm device handle is not rotated by ensuring bubble of level is between black lines (**Figure 7**) and insertion tube is not flexed throughout Blended Polymer delivery.
- d. With same free hand, push plunger forward **SLOWLY** (over 5 seconds) until it clicks and stops (approximately half plunger length) indicating both balloons are fully inflated (**Figure 9**). Hold device steady for 15 seconds.

Troubleshooting

Plunger will not advance during balloon inflation

Ensure slider is completely forward indicating catheters are fully advanced.

Figure 5


Figure 6


Figure 7


Figure 8



Figure 9

FDS-200 • Delivery System

FemBloc® Permanent Birth Control

INSTRUCTIONS FOR USE

- e. With same free hand, slide lock to the right until it clicks (*Figure 10*).
- f. With same free hand, push plunger forward **SLOWLY** (1 ratchet click per 1 second) until plunger ring is completely flush with the handle indicating the entire volume of Blended Polymer has been delivered (*Figure 11*).

Note: Last few ratchet clicks may require more force to advance plunger; ensure steady grip on device handle to prevent forward movement while advancing plunger.

Note: Avoid intermittent delivery of Blended Polymer, which may result in Blended Polymer solidifying at tip of catheter(s), preventing complete delivery.

- g. Hold device steady for 15 seconds.

CAUTION: Holding Delivery System in patient for longer than 15 seconds may lead to patient discomfort or injury during device removal.

- h. Detach tenaculum stabilizer while holding device handle steady.

Figure 10


Figure 11


Figure 12


Figure 13



Figure 14

IV. COMPLETE PROCEDURE

1. REMOVE DELIVERY SYSTEM

- a. With same free hand, slide lock to left until it clicks (*Figure 12*).
- b. Slowly retract plunger with free hand to deflate both catheter balloons (approximately half plunger length) (*Figure 13*).
- c. Move slider back to retract catheters (*Figure 14*).
- d. Remove device and discard.

Note: Solidified Blended Polymer may appear on insertion tube and/or balloons.

- e. Detach tenaculum and control for any bleeding at tenaculum site.

2. POST PROCEDURE INSTRUCTIONS AND PATIENT FOLLOW-UP

- a. Ensure patient uses alternative contraception (except IUD) until the Confirmation Test. Counsel patient on risk of pregnancy due to non-compliance.
- b. Schedule patient for Confirmation Test at least three (3) months after the FemBloc procedure.

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
Tel: (31) (0) 70 345-8570

Manufacturer

Femasys Inc.

3950 Johns Creek Court, Suite 100

Suwanee, GA 30024 USA

Tel: (877) 336-2562

Fax: (404) 581-5903

www.femasys.com

Patent Listing: www.femasys.com/patents

www.fembloc.com

VIGILANCE REPORTING

Serious incidents occurring in the European Union in relation to FemBloc Delivery System should be reported to the Authorized Representative and/or the Manufacturer as indicated below:

Authorized Representative: Emergo – EmergoVigilance@ul.com

Manufacturer: Femasys Inc. – Vigilance@femasys.com


User/Patient: Serious incidents should be reported to the manufacturer and corresponding competent authority of the Member State.

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Fabricante: Femasys Inc. - Vigilance@femasys.com